



# ENGINEERED FOR A NEW STANDARD. YOURS.

You know what the right outcomes look like.  
Rapid first-pass intubation. Lower risk of complications.  
Improved patient care. Introducing the video laryngoscopes  
designed to make them a reality—for both routine  
and difficult airways.



# GLIDESCOPE SINGLE-USE

## VIDEO LARYNGOSCOPES



# READY. JUST LIKE YOU.

When it comes to patient care, you don't leave anything to chance. Neither do we. That's why new GlideScope Single-Use Video Laryngoscopes have been so rigorously designed for your world. With a new ergonomic design and the industry's slimmest-profile blades, GlideScope gives you the maneuverability and working space you need to intubate adult patients quickly and safely, in routine and difficult airways.

## Low profile. High performance. The new LoPro S3 and S4 blades

### Focus on your patients instead of your equipment.

Our lowest-profile single-use blades—now featuring an ergonomic handle design for added control and comfort—give you more maneuverability and working space in the airway during even the most difficult intubations.



**ERGONOMIC HANDLE**

LoPro S3  
height at camera:  
**11.0 mm**

**LOW PROFILE DESIGN**

for improved maneuverability  
and working space

**SIGNATURE GLIDESCOPE  
BLADE ANGULATION**

improves your view of the glottis

**DIGITAL CAMERA  
LOCATION**

provides a wider forward view



### GlideScope Video Monitor

Get a better airway view with a large, full-color video monitor. Confirm tube placement and improve training with features such as snapshot, video recording and onboard tutorial.



### GlideScope Smart Cable

Discover a more elegant and versatile clinical solution with our exclusive GlideScope Smart Cable, which houses miniaturized video electronics in a titanium casing, that converts images to video. Durable and easy to clean and disinfect. The Smart Cable is the foundation of the platform and can be used with new and future single-use products.

## The single-use advantage

Get superior convenience without sacrificing the advantages of reusable blades. Sterile-packed and ready-to-use, GlideScope Single-Use blades help reduce infection opportunity and limit the need for instrument cleaning and reprocessing.

## Comfort, meet confidence The new MAC S3 and S4 blades


### The familiarity of Macintosh-style blades meets the performance of GlideScope.

Our new GlideScope Single-Use MAC blades give you the feel of the traditional handle you're used to working with. And with the superior view into airways that video laryngoscopy provides, you'll know more about your patient's airways—and can act on that existing knowledge.

**RECESSED CAMERA**  
helps prevent blood and secretions  
from obstructing your view

**ANTI-FOG**  
for a clearer view

MAC S3  
height at camera:  
**11.2 mm**




**TRADITIONAL  
MAC-STYLE BLADE**


## THE RIGHT OUTCOMES

In a recent study of novice users, video laryngoscopy improved success rate and time-to-intubate when compared with direct laryngoscopy with a MAC blade.<sup>1</sup>

### FIRST-PASS SUCCESS:



### TIME-TO-INTUBATE:




1. Ambrosio, Art, et al. "Difficult Airway Management for Novice Physicians: A Randomized Trial Comparing Direct and Video-Assisted Laryngoscopy." Otolaryngology—Head and Neck Surgery 150.5 (2014): 775-778.

# THE TOOLS YOU TRUST. FROM THE BRAND YOU KNOW.

Every GlideScope instrument we put in your hands is backed by years of industry leadership and engineered for the real situations you face when caring for patients. From new blade designs to the precision technology today's clinical settings demand, Verathon provides real-world innovation for reliable outcomes.

Ask your Verathon representative about **GlideScope Titanium Single-Use**, and explore our full range of products at [verathon.com](http://verathon.com).

## The GlideScope Video Laryngoscope Family



Corporate Headquarters  
20001 North Creek Parkway / Bothell, WA 98011  
425.867.1348 / 800.331.2313 / Fax: 425.883.2896  
[verathon.com](http://verathon.com)

**EC REP** Verathon Medical (Europe) B.V.  
Linnaeusweg 11 3401 MS IJsselstein The Netherlands  
+31.30.68.70.570 / Fax: +31.30.68.70.512  
[verathon.eu](http://verathon.eu)

GlideScope, the GlideScope symbol, Verathon, and the Verathon Torch symbol are trademarks of Verathon Inc. © 2014 Verathon Inc.

CE 0413

0900-4579-00-86